

# photomultiplier bench top HV supply

## HVLAB3000 data sheet

### 1 description

The HVLAB3000 is a programmable 3000V high stability, dc power supply for laboratory and test applications, designed with the photomultiplier user in mind. It has a small foot-print and features a choice of positive or negative polarity. Additional +/-5V and +/-12V auxiliary outputs are provided for powering associated electronics modules such as amplifiers and photon counting instrumentation. Power input is 100V to 240V ac.

The output voltage can be set using three different methods: manually, using a ten-turn precision front-mounted potentiometer, digitally via a USB port, or using an analogue control input. Two power supplies can be controlled by the HVLAB3000 software.

The voltage setting is retained on switch-off and upper limits for HV and current can be set manually via the rear panel adjustment potentiometers using a small screwdriver (tamper proof). The current limit folds back the high voltage to limit the current. The voltage limit is the same for both polarities and the current limits are set individually for each polarity. Both voltage and current limits operate independently of the soft-set control limits.

Adjustable brightness LED displays provide accurate readings of the HV output voltage and current. The same displays indicate the set maxima for HV and current by push button, as selected by the user.

The USB control is compatible with the ET Enterprises MCS-CT3 multi-channel scaler when used in photomultiplier photon counting systems.

### 2 applications

- general laboratory use
- photomultiplier based photon detection systems

### 3 features

- compact
- HV adjustable from 100-3000V
- user selectable HV polarity
- overload protected
- adjustable brightness LED displays
- manual or USB controlled HV adjustment
- user selectable upper current limit
- small footprint
- dual polarity auxiliary low voltage output
- USB interface
- analogue remote control (0-3V)






### 4 specifications

	unit	min	typ	max
<b>input voltage (50/60Hz)</b>	V	100		240
<b>remote control input</b>	V	0		3
<b>input power (full load)</b>	VA		46	
<b>output voltage (+ or -)</b>	V	100		3000
<b>output ripple (full load)</b>	% full scale		0.01	
<b>output current</b>	mA			3
<b>voltage set precision:</b>				
manual	V		1	
software	V		1	
<b>voltage regulation</b>	%			0.001
<b>load regulation (0-2.5mA)</b>	%		1	
<b>discharge time (90 to 10% with 4M7 load)</b>	s		4	
<b>auxiliary voltage output 1:</b>				
voltage (fixed)	V		+/-5	
current	mA			300
<b>auxiliary voltage output 2:</b>				
voltage (fixed)	V		+/-12	
current	mA			125
<b>protection</b>				
short circuit (to ground)		protected against continuous short circuits		
<b>temperature</b>				
operating	°C	0		40
storage	°C	-40		60
weight	kg		1.5	

## 5 connections

input connector	IEC SHV
HV output connector (2)	SHV
interface connector	USB B (2.0)
auxiliary output connector	DIN
remote control sockets (2)	4mm banana socket

## 6 outline dimensions (mm)



## 7 ordering information

item	ordering code
3000V dual polarity HV cable to PMT housing	HVLAB3000 LEADMHVSHV-100cm

## 8 warning

High voltages generated by these products present an electrical shock hazard and appropriate precautions must be taken. Installation must be by qualified personnel.

All units are despatched with the HV control potentiometer set to zero.

Do not operate outside the specification limits of the HVLAB3000 or those of the photomultiplier. This may result in loss of performance, permanent damage, or both.

**ET Enterprises Limited**  
45 Riverside Way  
Uxbridge UB8 2YF  
United Kingdom  
tel: +44 (0) 1895 200880  
fax: +44 (0) 1895 270873  
e-mail: sales@et-enterprises.com  
web site: www.et-enterprises.com

**ADIT Electron Tubes**  
300 Crane Street  
Sweetwater TX 79556 USA  
tel: (325) 235 1418  
toll free: (800) 399 4557  
fax: (325) 235 2872  
e-mail: sales@electron tubes.com  
web site: www.electrontubes.com

an ISO 9001 and ISO 14001 registered company

The company reserves the right to modify these designs and specifications without notice. Developmental devices are intended for evaluation and no obligation is assumed for future manufacture. While every effort is made to ensure accuracy of published information the company cannot be held responsible for errors or consequences arising therefrom.

ET Enterprises  
electron tubes

© ET Enterprises Ltd, 2016  
DS\_ HVLAB3000 Issue 2 (11/02/16)